
Spare Part Replacement Manual for **Ambu[®] aView[™] 2 Advance**

Content

Page

English	3
България	31
Česky	60
Dansk	89
Deutsch	118
Ελληνικά	147
Español	176
Eesti	205
Suomi	234
Français	263
Hrvatski	292
Magyar	321
Italiano	350
日本語	379
Lietuviškai	408
Latviski	437
Nederlands	466
Norsk	495
Polski	524
Português	553
Romania	582
Русский	611
Slovenčina	640
Slovenščina	669
Svenska	698
Türkçe	727
中文	756

Ambu is a registered trademark and aScope and aView are trademarks of Ambu A/S.

Content

Page

1. Important information	4
1.1. Requirements and Qualifications for Personnel Performing Spare Part Replacement	4
1.2. Warnings and Cautions	4
1.3. Data security	4
1.4. System Description	5
1.5. Spare Parts	6
1.6. System Specifications	6
1.7. Spare Part Replacement Process	8
2. How To Clean and Disinfect the Displaying Unit	9
3. Preparing for Spare Part Replacement	10
3.1. Tools Needed	10
3.2. Work Environment	10
3.3. How to Disassemble the Displaying Unit	10
4. How to Replace the Spare Parts	15
4.1. Battery Replacement	15
4.2. Visualization Device Interface Board Replacement	17
4.3. Fan Replacement	21
4.4. Insert a New USB Cover	23
4.5. How To Reassemble the Displaying Unit	24
5. Testing	28
5.1. Visual Inspection	28
5.2. Protective Earth Resistance Test	28
5.3. Insulation Resistance Test	28
5.4. Leakage Current Tests	29
5.5. Functional Tests	29

1. Important information

This is a Spare Part Replacement Manual for Ambu® aView™ 2 Advance. It is applicable to aView 2 Advance with item number 405011000 ONLY.

The following terms are used in this Spare Part Replacement Manual:

- “Displaying unit” refers to Ambu aView 2 Advance.
- “Visualization device” refers to the various Ambu scopes that can be connected and used with the displaying unit.

This manual may be updated without further notice. Contact your local Ambu representative or download the most recent version from www.ambu.com.

Please observe that the warranty will be void, if the displaying unit has been disassembled within the warranty period without the written consent from Ambu.

1.1. Requirements and Qualifications for Personnel Performing Spare Part Replacement

Do not start replacing any spare part of the displaying unit, unless you have read through this document.

Do not start replacing any spare part of the displaying unit, unless you are qualified in the following areas:

- Knowledge, experience and acquaintance with electronic repair and safety testing according to IEC 60601-1 (Medical electrical equipment, Part 1: General requirements for basic safety and essential performance) and IEC 62353 (Medical electrical equipment – Recurrent test and test after repair of medical equipment).
- Knowledge of and experience with local regulations.
- Having read the Instructions for Use (IFU) for the displaying unit (download current version from www.ambu.com).
- Knowledge of the environment where the displaying unit is installed and used.

1.2. Warnings and Cautions

WARNINGS

- To avoid risk of electric shock, always disconnect the mains from the displaying unit before performing spare part replacement.
- To avoid contamination, always clean and disinfect the displaying unit before and after the replacement and wear gloves during handling.
- To avoid compromising electrical safety, make sure not to leave any tools or other loose parts inside the displaying unit when it is reassembled.
- To avoid compromising electrical safety, always perform tests after the displaying unit has been reassembled. If any part of the tests fails after spare part replacement, do not use the displaying unit.
- Only use spare parts supplied by Ambu. Do not modify the spare parts.
- For battery spare part kit: Risk of fire and burns. Do not open, crush, heat above 45 °C or incinerate the battery.

CAUTIONS

- Use ESD protection when disassembling, replacing any spare part, or reassembling the displaying unit.
- The displaying unit must be turned OFF before disassembly.
- Before starting the replacement of any spare part, make sure that the battery is disconnected.

1.3. Data security

To avoid possible loss of data, it is recommended to export relevant files before performing any spare part replacement. Refer to the online Reference Manual at www.ambu.com.

1.4. System Description

The displaying unit is a non-sterile, reusable digital displaying unit, intended to display live imaging data from Ambu visualization devices.


Figure 1 System illustration

No.	Part	Function
1	Displaying unit	Display live imaging data from Ambu visualization devices.
2	Visualization device interface (VDI) board	The interface board that connects to Ambu® visualization devices.
3	LCD Touch Screen	Displays the image from the Ambu visualization device and provides users with graphical interface.
4	Colour Ring	Indicates which type of visualization device it supports. The colour of the ring must match the colour on the connector of the visualization device.
5	Visualization device connector cover	Currently not in use.
6	Stand	Use the stand to place the displaying unit on a solid surface or to carry the displaying unit while turned OFF.
7	Power button	Push button for power ON and OFF.
8	Hardware reset button	Resets the displaying unit hardware without impacting stored data.
9	Fan	Cooling the system.
10	Power inlet	Power inlet for charging the displaying unit.
11	Input/output connections	HDMI, SDI, LAN.
12	2 x USBs port with USB cover	Connects to a USB memory device for file export or software update.
13	Battery	Powers the system.
14	Bracket	Secures the displaying unit to e.g. an IV pole.
15	Power supply bracket	Secures the placement of the power supply.





No.	Part	Function
16	Release buttons	Releases displaying unit from the bracket.
17	Hook	Holder for pouches containing visualization devices.
18	Power supply	Powers the displaying unit.
19	Power cable	Power cable with country-specific plug.

1.5. Spare Parts

The spare parts may not be available in all countries. Please contact your local Ambu representative.

WARNING

Only use spare parts supplied by Ambu. Do not modify the spare parts.

Spare Part	Description	Item numbers:
	Ambu® aView™ 2 Advance - Bracket	405013100
	Ambu® aView™ 2 Advance - Power Supply Kits: I (AUS, ZH) B (US, JP) E / F / J (EU – Not for DK, UK) G (UK) K (DK)	405015100 405015200 405015300 405015400 405015500
	Ambu® aView™ 2 Advance - Battery Kit This kit contains a battery and a replacement part kit. See section 4.1.	405012100
	Ambu® aView™ 2 Advance - Visualization Device Interface Kit - Blue This kit contains a visualization device interface board, a colour ring and a replacement part kit. See section 4.2.	405014100
	Ambu® aView™ 2 Advance - Fan Kit. This kit contains a fan and a replacement part kit. See section 4.3.	405016100

1.6. System Specifications

Specifications for aView 2 Advance

Display	
Resolution	1920 x 1080 pixels
Orientation	Landscape
Display type	12.8" colour TFT LCD

Connections	
2 x USB connection	USB 3.0 Type A
Digital video out	HDMI and 3G-SDI (1920 x 1080 p, 60 fps)
WIFI	Supports IEEE Standards 802.11ac/a/b/g/n
LAN	Ethernet RJ45 connector, 10/100/1000 Mbps
Memory	
Storage capacity	32GB
Mounting interface	
Mounting interface standard	VESA MIS-D, 75 C, VESA FDMI compliant display, Part D, with centre located mounting interface
Electrical power	
Power requirement	19.0 VDC; 3.43 ADC
Battery type	14.4 VDC 6500mAh
Battery Operation	Typical battery runtime of a new, fully charged battery (the displaying unit turned on and scope connected) is more than 3 hours
Protection against electrical shock	Internally Powered
Environment conditions	
Transportation temperature	Minus 10 - 55°C (14 - 131°F)
Storage temperature	10 - 40°C (50 - 104°F)
Operation temperature	10 - 35°C (50 - 95°F)
Relative humidity	30 - 85%
Atmospheric pressure	80-109 kPa
Altitude	≤ 2000 m
IP Protection Classification System	The aView 2 Advance is classified IP30. Protection against solid objects
Dimensions	
Width	331 mm (13.03")
Height	215 mm (8.46")
Thickness	52 mm (2.05")
Weight	2.7 kg (6.0 lbs)
Bracket	
Fits poles with thicknesses	Ø 18 mm - 35 mm (0.7 - 1.5")

Specifications for aView 2 Advance Power Supply

Dimensions	
Weight	360g (0.79 lbs)
Electrical power	
Power requirement	100 - 240V AC; 50-60 Hz; 1.0-2.0 A
Power out	19.0VDC 3.43ADC
Protection against electrical shock	Class I

Operating and storage environment	
Temperature	10 - 40° C (50 - 104° F)
Relative humidity	30 - 85%
Plugs	
Between the power supply and the displaying unit	Ø6.5mm DC jack connector
6 interchangeable types	USA and Japan: Model NEMA 5 AC grounded power plug Australian: AS3112, AC grounded power plug UK: BS1363, AC grounded power plug European: CEE 7, AC grounded power plug Danish: 2-5a, AC grounded power plug Swiss: Type J, AC grounded power plug

1.7. Spare Part Replacement Process

This section specifies a spare part replacement process, which is in compliance with IEC 62353.

The organisation and personnel being responsible for the process should fulfil the requirements outlined in section 1.1. Follow the process in Figure 2. to replace spare parts in the displaying unit.


Figure 2 Spare part replacement flow

2. How To Clean and Disinfect the Displaying Unit

The displaying unit shall be cleaned and disinfected before and after performing replacement according to good medical practice, described in the instructions below.

The displaying unit can be cleaned and disinfected using disposable germicidal hospital disinfection wipes designed for hard, non-porous surfaces. Compatible disinfection wipes are based on Isopropyl alcohol 70-80%, Isopropyl alcohol/Ammonium Chloride or Isopropyl alcohol/Ethanol.

- a) Turn OFF the displaying unit and make sure that it is unplugged from the power outlet.
- b) Follow the cleaning instruction from the germicidal disinfection wipe manufacturer to thoroughly clean the displaying unit.
- c) Allow the displaying unit to dry after the cleaning.
- d) Follow the disinfection instruction from the germicidal disinfection wipe manufacturer to thoroughly disinfect the displaying unit with a new wipe.
- e) Allow the displaying unit to dry after the disinfection.
- f) If the cleaning and disinfection chemicals used leave a residual layer on the display after they dry out, use an alcohol-based wipe to remove the residue.
- g) If applicable, clean and disinfect the mounting bracket following the same procedure as for the displaying unit.

NOTES:

- The wipes shall be moist, but not dripping to ensure no damage to internal electronics of the displaying unit.
- If using wipes containing hypochlorite, ensure that all residue is completely removed. Wipes containing hypochlorite may affect the screen's antireflective coatings over time. You should limit the use of wipes containing hypochlorite to required cases only.

3. Preparing for Spare Part Replacement

3.1. Tools Needed

- Name plate removal tool
- Plectrum tool
- Screwdriver(Phillips #0) – Not supplied
- Compatible (and functioning) Visualization Device – Not supplied

The name plate removal tool and the plectrum tool are provided with the spare parts kits (see section 1.5 Spare Parts).


Figure 3 1 Name plate removal tool, 2 Plectrum tool

3.2. Work Environment

CAUTION

Use ESD protection when disassembling, replacing any spare part, or reassembling the displaying unit.

It is recommended to work on a soft surface to avoid scratches on the screen and surface of the unit.

3.3. How to Disassemble the Displaying Unit

WARNING

To avoid risk of electric shock, always disconnect the mains from the displaying unit before performing spare part replacement.

CAUTION

The displaying unit must be turned OFF before disassembly.

Follow these steps to disassemble the displaying unit:

a) Clean the displaying unit before disassembly

WARNING

To avoid contamination, always clean and disinfect the displaying unit before and after the replacement and wear gloves during handling. Follow the instructions provided in section 2 How To Clean and Disinfect the Displaying Unit.

b) Remove all cables attached to the displaying unit

E.g. power supply, visualization devices and external equipment

c) Remove the name plates (with the Ambu logo) on the top and bottom of the displaying unit

Refer to Figure 4: Insert the name plate removal tool in the hole at the end of the name plate. Press the name plate removal tool down to flip the end of the name plate up. Use the tool as a lever. Gently lift the name plates upwards.


Figure 4 Removing the name plate

d) Remove the round screw covers on each side of the displaying unit

Insert the name plate removal tool in the hole by the screw cover. Press the name plate removal tool down to lift the screw cover.

Figure 5 Removing the screw cover

e) Remove the screws on all four edges of the displaying unit

Use the screwdriver to remove the screws. Total 10 screws must be removed: 4 screws on top, 4 on the bottom, and 1 screw on each side.


Figure 6 Removing the screws

f) Remove the back casing

Use the plectrum tool. Press it into the dividing line between the displaying unit and the back casing. On all four edges, start in the middle and pull the plectrum tool towards the corners until the back casing is released. When the two sides of the casing is disassembles, you will here a “click” sound, indicating that the locks are released. (Figure 7)

Figure 7 Releasing the back casing using the plectrum tool

- g) Place the displaying unit on a flat and soft surface facing down (see Figure 8).**

Figure 8 Place the displaying unit facing down

- h) Lift the back casing off**

At the opposite side of the colour rings, as indicated in Figure 8, carefully pull the back casing 2 cm upwards (see Figure 9). Do not move the back casing too far before the internal FPC cable is released. Slide the back casing sideways 1 cm (see Figure 9). Do not move the back casing too far, as it may accidentally damage the internal FPC cable.


Figure 9 Carefully pull the back casing upwards and slide it to the right slightly

Carefully open the back casing, while the internal FPC cable is still connected.

Figure 10 Open the back casing carefully and remove adhesive tape where necessary

i) Remove adhesive tape

If adhesive tape is protecting the connectors on the boards, you can remove it (see Figure 10). It is not necessary to replace the tape, as it is only for protection during transportation to the end user.

j) Release the FPC cable

Flip the white lock upwards with a fingernail (see Figure 11).


Figure 11 Disconnect the FPC cable

With the FPC cable released, lift the back casing and put it on the table.
Now the main board is exposed, see Figure 12.


Figure 12 Overview of the main board. Part names: 1 Fan connector, 2 FPC connector to main board, 3 Fan, 4 Battery connector, 5 Battery, 6 Visualization device Interface (VDI) Board 1, 7 Visualization device Interface (VDI) Board 2.

k) Disconnect the battery cable

CAUTION

Before starting the replacement of any spare part, make sure that the battery is disconnected.

Release the battery cable (see location on Figure 12) by pulling the grey part carefully away from the black part in the direction shown on Figure 13.


Figure 13 Disconnect the battery cable

4. How to Replace the Spare Parts



This section covers the replacement of different parts. You may jump directly to the section covering the specific part to be replaced:

- Battery (section 4.1)
- Visualization device interface (VDI) board (section 4.2)
- Fan (section 4.3)
- USB cover (section 4.4)

WARNING

Only use spare parts supplied by Ambu. Do not modify the spare parts.

4.1. Battery Replacement

Ambu® aView™ 2 Advance - Battery Kit	405012100
Battery	
Replacement part kit: <ol style="list-style-type: none"> 1 Name plates (x2) 2 Screw cover (x2) 3 Visualization device connector cover (x2) 4 Plectrum tool (x1) 5 USB Cover (x1) 6 Screws for casing assembly (x10) 7 Screws for VESA the mounting bracket (x4) 8 Name plate removal tool (x1) 	

WARNING

Risk of fire and burns. Do not open, crush, heat above 45 °C or incinerate the battery.

Replacing Battery

Refer to Figure 12 to identify the location of the battery.

- a) **Pull the battery off the displaying unit main board**
The battery is fastened with double-sided adhesive tape.

Figure 14 Remove the battery

b) Dispose the used battery according to local guidelines for disposal of Li-ion batteries

c) Remove the remaining tape from the displaying unit

Remove as much as possible of the old double-sided adhesive tape.

d) Place the new battery




Remove the cover on the double-sided tape on the back of the new battery. The battery must be placed with the labels pointing to the right and on the upper side of the battery (see Figure 15).

Figure 15 Place a new battery, with the adhesive tape facing down

e) Follow the steps in chapter 4.5 Reassembly of the displaying unit

4.2. Visualization Device Interface Board Replacement

Ambu® aView™ 2 Visualization Device Interface Kit - Blue	405014100
VDI board - Blue	
Colour ring	
Screws for VDI board (x2)	
Replacement part kit: 1 Name plates (x2) 2 Screw cover (x2) 3 Visualization device connector cover (x2) 4 Plectrum tool (x1) 5 USB Cover (x1) 6 Screws for casing assembly (x10) 7 Screws for VESA the mounting bracket (x4) 8 Name plate removal tool (x1)	

Additional Tools (not supplied)

Screwdriver (Phillips #0)

Replacing the VDI board

Refer to Figure 12 to identify the location of the VDI board.

- a) Use the screwdriver to loosen the two black screws holding the VDI board (see Figure 16)


Figure 16 Remove the screws on the VDI board to be replaced

- b) Pull the VDI board out of the main board connector


Figure 17 Remove the VDI board

- c) **Dispose the used VDI board following local guidelines for disposal of electronic waste**
- d) **Insert the new VDI board into the main board connector**


Figure 18 Insert a new VDI board

- e) **Insert and tighten the two black screws (max. torque 0.2Nm) holding the VDI board in place (see Figure 16)**
- f) **Make sure the colour ring is the same colour as the indication on the VDI board (see Figure 19)**
When remounting the back casing, be aware that the coloured spots on the VDI boards (Figure 19) are matching the colours of the colour rings at which they are placed.


Figure 19 The colour indication on the VDI board must match the colour ring

Replace a damaged colour ring

If the colour ring is broken, replace it with the new colour ring that comes with the VDI kit.

- a) **Remove the damaged colour ring by pushing it out.**


Figure 20 Remove the damaged colour ring

- b) **Push the new colour ring into place from the outside of the back casing**

Make sure that the squared cut-out in the colour ring is fitted to the corresponding shape on the inside of the back casing (see Figure 21).


Figure 21 Insert a new colour ring

Inserting a visualization device connector cover, if necessary

If a VDI board has been removed from its position, it is necessary to cover the hole in the casing with a visualization device connector cover.

a) Insert a visualization device connector cover



Push the visualization device connector cover in place from the outside of the back casing (see Figure 22).

Figure 22 Install the cover if, for any reason, the VDI board has been removed

b) Follow the steps in chapter 4.5 to reassemble the displaying unit

4.3. Fan Replacement

Ambu® aView™ 2 – Fan Kit	405016100
Fan	
Screws for fan Replacement part kit (x3)	
Replacement part kit: 1 Name plates (x2) 2 Screw cover (x2) 3 Visualization device connector cover (x2) 4 Plectrum tool (x1) 5 USB Cover (x1) 6 Screws for casing assembly (x10) 7 Screws for VESA the mounting bracket (x4) 8 Name plate removal tool (x1)	

Additional Tools

Screwdriver (Phillips #0)

Replacing fan

Refer to Figure 12 to identify the location of the fan.

a) Disconnect the fan connector

Identify the location of the fan connector on the main board (see Figure 12), and disconnect the fan connector (see Figure 23) .


Figure 23 Disconnect the fan cable

b) Untighten the three screws holding the fan in place (see Figure 24)


Figure 24 Remove the screws from the fan

- c) **Dispose the used fan according to local guidelines for disposal of electrical waste**
- d) **Insert the new fan**
- e) **Install the new fan with the cable pointing towards the upper edge of the displaying unit.**

Figure 25 Install the fan

- f) **Insert and tighten the three screws (max. torque 0.3Nm) holding the fan in place (see Figure 24)**
- g) **Connect the fan cable (see Figure 12 and Figure 23)**
- h) **Follow the steps in chapter 4.5 to reassemble the displaying unit**

4.4. Insert a New USB Cover

If the USB cover is damaged, a replacement cover is included in any of the spare part kits.

- a) Find and remove all parts of the the damaged USB cover from the inside of the displaying unit.
- b) Push the new USB cover in from the outside of the casing. You can using a flat screwdriver to assist but take care not to damage the rubber.


Figure 26 Install a new USB cover

- c) Follow the steps in chapter 4.5 to reassemble the displaying unit

4.5. How To Reassemble the Displaying Unit

WARNING

To avoid compromising electrical safety, make sure not to leave any tools or other loose parts inside the displaying unit when it is reassembled.

a) Reconnect the battery to the battery connector (see Figure 12 and Figure 27)

Ensure that the connector is inserted properly by observing that the white centre guide is touching the black connector part


Figure 27 Reconnect the battery connector

b) Reconnect the FPC cable

Reconnect the FPC cable to the main board. Also ensure the cable is still connected to I/O board on the back casing. The following figure shows the location of FPC cable and FPC connectors.


Figure 28 The locations for FPC cable and FPC connectors.

- 1 FPC (Flexible Printed Circuit) cable,
- 2 FPC connection to the main board
- 3 FPC connection to I/O board on the back casing.

Reconnect the FPC cable to the main board (see location 2 in Figure 28). Insert the FPC cable in the connector and lock the white flat by flipping it down (see Figure 29). Make sure the golden leads on the printed circuit are facing down and have good contact with the connector.

Make sure that the end of FPC cable is pushed all the way and parallel to the connectors. Check that the FPC cable is still connected correlly to the I/O board (See location 3 in Figure 28).


Figure 29 Connect the FPC cable back to the main board

- c) **Put on the back casing, starting with the side where the colour rings are placed (see Figure 30).** Make sure that the edge of the back casing is aligned and pushed against the edge of the displaying unit front underneath the colour rings.

Figure 30 Place the back casing, starting from the side with the colour rings

- d) **Pull the back casing towards and around the USB ports and push the back casing in place (see Figure 31 and Figure 32).**

Make sure that the edge of the back casing does not damage the USB ports. Make sure that all the locks around the edge of the displaying unit are clicked in place.

Figure 31 Put on the back casing, continued

Figure 32 Put on the back casing till it clicks in place

e) Re-insert screws on all four edges of the displaying unit

Insert 10 screws in total around the displaying unit edge. Four screws on top and bottom, and one screw on each side. If any screws are missing use the new ones supplied with the spare part kit.

f) Click the name plates in place (see Figure 33)

Click down the two ends of the name plate first (see Figure 33). Then click down the rest of the name plate towards the center. The Ambu logo on the name plate must be readable from the front side of the displaying unit. If the name plates are damaged, use the new ones supplied with the spare part kit.


Figure 33 Put on the name plates

g) Click the screw covers in place (see Figure 34)

When clicking the screw covers in place, be aware that they must be oriented correctly, in order to fit (see Figure 34). If the screw covers are damaged, use the new ones supplied with the spare part kit.

Figure 34 Put on the screw covers

h) Perform safety test

Follow the instructions provided in 5 Testing.

WARNING

To avoid compromising electrical safety, make sure not to leave any tools or other loose parts inside the displaying unit when it is reassembled.

i) Clean the displaying unit after reassembly and testing

Follow the instructions provided in section 2 How To Clean and Disinfect the Displaying Unit.

5. Testing

The testing instructions in this chapter comply with the requirements of IEC 62353.

WARNING

To avoid compromising electrical safety, always perform tests after the displaying unit has been reassembled. If any part of the tests fails after the spare part replacement, do not use the displaying unit.

When the displaying unit has been disassembled, it is mandatory to perform all the tests specified in the following sections.

5.1. Visual Inspection

During the process of replacing the spare parts, the inner parts shall be inspected for any damage, contamination or spilling.

When the displaying unit has been re-assembled, make sure that the ID marking on the back of the unit is complete and clearly visible, and that all warnings and marking icons by the ports are complete and clearly visible.

Inspect if any damage, contamination or spilling may have compromised the safety of the displaying unit.

Inspect if all mechanical parts are in good shape and not broken.

Inspect if the power supply and its cords are in good shape and not broken.

5.2. Protective Earth Resistance Test

Not applicable

As the displaying unit is not directly connected to the supply mains (110/230 V AC), the Protective Earth Resistance test is only relevant for the power supply unit.

Since the Power Supply unit is not a repairable part, this test is not applicable to the displaying unit.

5.3. Insulation Resistance Test

Insulation resistance test according to 62353 is required after replacing the spare parts of the displaying unit.

5.4. Leakage Current Tests

Leakage current test is required after replacing the spare parts. When testing Applied Part Leakage Current, an Ambu visualization device shall be inserted, since the visualization device is the applied part.

In order to be able to measure applied part leakage current, use metal foil wrapped around the distal tip of the visualization device (See Figure 35).


Figure 35 Setup for the leakage current test

5.5. Functional Tests

Run any of the following functional tests relevant for the spare part replaced.

Image Display and Visualization Device Connector Test

After the displaying unit has been re-assembled, always perform a test of the image to check if the replaced parts and the whole system works as expected.

- a) Turn the displaying unit ON
- b) Connect a visualization device and point the distal end of the visualization device towards an object, e.g. the palm of your hand.
- c) Verify that a live video image appears on the screen
- d) Make sure that the live image has the correct orientation
- e) Check that there are no lines or other unexpected disruptions in the image
- f) Check that there is light in the visualization device LED
- g) Repeat steps II-VI for each of the visualization device connectors
- h) Verify the touch function by touching some of the options on the screen.

Battery Check

Check that the battery discharges as expected.

- a) Connect the power supply and charge the displaying unit until the battery is fully charged
- b) Disconnect the power supply
- c) Wait for the screen saver to be activated (as this will change the expected remaining battery time)
- d) After 10 minutes check the time in the battery icon.

If the battery remaining time decreases as expected $\pm 20\%$, the test is passed (e.g. battery remaining time goes from 4:10 to anywhere between 3:58 to 4:02 in 10 min.).


Fan Check

- a) Turn off the displaying unit.
- b) When turning on the unit, check if the fan starts rotating silently by looking through the ventilation ribs and listening for any noise.

External Monitor Connector Check

This test checks if the FPC cable (see Figure 12) has been connected correctly, so the external monitors can be used.

- a) Plug in an external monitor in the HDMI or the SDI plug.
- b) Check if the image is shown on the monitor in good quality.

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